


General Comments

Dunlop Foams recommend the mattresses are used on beds which have mattress retainers where users are concerned about the possibility of a mattress sliding off the surface of a particular bed.



The distance between the ends of the mattress and the head / foot boards should never be more than 30mm.

Non perforated mattress platforms on electric beds may inhibit air circulation around the mattress base. This may result in the mattress maintaining heat for a longer period. Please note that patients in an environment which is not climate controlled / air conditioned may sweat, which left unchecked, could contribute to the development of a pressure injury.